

Accessories & Spares

Airway, Guedal, Disposable	Order No
Paediatric Size 00	W6873
Paediatric Size 0	W6874
Child Size 1	W6875
Adult Size 2	W6876
Adult Size 3	W6877

Carrying Options	Order No
Bags - see page 36 - 39 for options	
Knapsack	W8654
Pouch, Belt Carrying	500-170

Control Module	Order No
Model 2R	500-A313/CE

Cylinders, Oxygen, Aluminium	Order No
Comprehensive Size 'D' 365L - Full	W6837
Comprehensive Size 'D' 365L - Empty	W6838
Compact Size 'C' 230L - Full	W6833
Compact Size 'C' 230L - Empty	W6834

Hoses, Input	Order No
Input Hose from Regulator to Control Module, 2 Metres long, for Model 2R, 4R and 4R(R) sets only, DISS to BS5682 O ₂	500-A226/2.00
Input Hoses up to maximum 15 metres available - add to above price	/Length

Hoses, Output	Order No
Output Hose from Control Module to Patient Valve, 1.5 metres long, for Model 2R, 4R & 4R (R) sets only	500-A200
Output Hoses up to maximum 2 metres available - add to above price	/Length

Hoses, Oxygen Therapy	Order No
Schrader BS5682 O ₂ - 6 Litre	500-A254/6

responder

Accessories & Spares Continued

Masks, Face	Order No
Adult Clear	W6807
Child/Small Adult Clear	W6806
Paediatric Size 0	W6803
Paediatric Size 1	W6804
Mouth Piece	W6808

Masks, Oxygen Therapy	Order No
Disposable Mask, Pack of 5	500-A369
Oxygen Tube (Pack of 5) for use with disposable therapy mask (A369)	500-A928/5

Miscellaneous Items	Order No
Exhaust Taper Fitting 30mm Male/22mm Female	500-A152
Gauge Contents	W6829CE
Harness, Clausen	500-437
Note with this option Hook Ring (W6872) will also be required	
Hook Ring	W6872
Labels, Location - Pack of 4	504-255
Lightstick Green	W6858
Lightstick Red	W6858/R
Nut, Mask Retaining - for 4R(R) only	500-145
Shroud, Gauge, Shock Absorbing	500-864
Spanner, Multi Purpose	500-A130
Washers, Pin Index Yoke	500-276
Wheel, Valve, Cylinder	500-A171

Oxygen Therapy Units	Order No
Schrader (incl. 5 disposable therapy face masks)	500-A253/6
Instantair (incl.5 disposable therapy face masks)	500-A263/3

Patient Pressure Monitor - Spares	Order No
Patient Pressure Monitor	500-A930/CE
Patient Valve Adaptor	500-A936
Threaded Adaptor	500-935
Ring & Tube Adaptor	500-A934
Connecting Hose (W6822 x 1.6m)	500-933
Filter (Pack of 5)	W8860/5
Gauge Housing	500-931/932
Gauge	W6831

responder

Accessories & Spares Continued

Regulators

Order No

See page 40 for options

Valves

Order No

Anti Inhalation Valve	500-A195/CE
O ₂ Therapy Control Valve (Schrader)	500-380/CE
Patient Valve without Relief Valve, Model 2R, 4R	500-A211/CE
Patient Valve without Relief Valve, Model 4R(R)	500-A212/CE
PEEP Valve, 4cm H ₂ O	500-A158/4
PEEP Valve, 8cmH ₂ O	500-A158/8
Relief Valve, Model 2R, 4R, 4R(R) with Audible Alarm, 40cm H ₂ O	500-A247/40/CE
Relief Valve, Model 2R, 4R, 4R(R) with Audible Alarm, 60cm H ₂ O	500-A247/60/CE

Users Manual

Order No

responder handbook	500-A251/CE
--------------------	-------------

Ordering Example: 500-A432CE/L6 + W1229 + W6838

- = **1)** Responder Model 2R (500-A432CE) with alternative Input Hose with AGA probe (suffix L6)
- + **2)** responder carrying bag (black/orange) (W1229)
- + **3)** 'D' size oxygen cylinder - empty (W6838)

Medical Deployment Bags

Instant Action Bags

Responder Bag

Easy to pick up and run with, the responder bag is a lightweight, purpose made carrying bag that can be hand carried, worn as a rucksack or over the shoulder. It has reflective strips for easy identification and neatly arranged compartments for a responder ventilator, cylinder, regulator, masks, Guedal airways and additional pockets for other medical supplies, providing instant access, whatever the circumstances.

PAC Bag

The PAC bag is a neat lightweight hospital carrying bag with reflective strips, hand, shoulder straps and metal hooks for hanging on a bed trolley. It contains a compartment to house a PAC ventilator with drive and patient hose, and comes with a replaceable base for maximum durability

PAC & Cylinder Bag

Ergonomically designed for paramedics the PAC & Cylinder bag provides instant access in emergency situations. It has compartments to house a PAC ventilator, drive hose, patient hose, cylinder, regulator and has additional pockets for other medical items such as laryngoscopes and suction units. This durable bag comes with a protective replaceable base, reflective strips and can be either hand carried, worn as a rucksack or over the shoulder.


Rucksack

This neatly designed backpack will carry any responder or PAC ventilator, regulator, cylinder, hoses and masks. It has specially designed pockets to house a defibrillator, suction, airway management and can be worn as a rucksack, leaving your hands free to carry any other emergency equipment. Reflective strips ensure easy identification.



Comprehensive Instant Action Bag

Designed for the Comprehensive Instant Action Set (see page 7), this red coloured hard wearing PVC coated nylon bag retains strength, durability and effectiveness while providing easy access, cleaning and exchange of cylinder. It is complimented by the day-yellow panels, making it easily visible and it has proven to provide for long service life. The case incorporates strong carrying handles and shoulder slings, and has Velcro fastenings for easy, quick and foolproof opening.



Compact Instant Action Bag

Smaller than the the comprehensive bag but equally effective this durable PVC coated nylon case is specially designed for the Compact Instant Action Set (see page 7). It incorporates strong carrying handles, shoulder sling and retains Velcro fastenings for rapid opening. The case can contain a responder or PAC ventilator, patient valve & lead, input lead, 230L O₂ cylinder and regulator.

Specification

responder Bag, PAC Bag, Oxygen & Cylinder Bag and Rucksack

Material: Heavy duty Cordura, washable and water resistant

Zips: Y2K high quality heavy duty zips

Labels: Removeable side labels with Velcro fastenings and clear pockets to slide in labelling. Flashed with Pneupac logo as standard. (Company) logo can be configured to customer requirements for which an additional cost may be charged.

Dimensions:

Responder Bag: 560mm W x 220mm H x 150mm D

PAC Bag: 360mm W x 200mm H x 160mm D

Oxygen & Cylinder Bag: 580mm W x 280mm H x 160mm D

Rucksack: 530mm W x 500mm H x 340mm D

Comprehensive Instant Action Case and Compact Instant Action Case

Material: Fire retardent PPO, with Velcro Fastenings

Dimensions:

Comprehensive: 625mm W x 300mm H x 300mm D

Compact: 515mm W x 150 mm H x 220 mmD

Ordering Information

Instant Action Bags	Order No
responder Bag	W1229
PAC Bag	W1231
PAC & Cylinder Bag	W1230
Rucksack	W1228
Comprehensive Instant Action Case (red)	500-A177
Compact Instant Action Case (grey)	500-A155

Medical Storage Bags



Trauma Pack Storage Bag

The Trauma Pack bag is a multiple storage bag, specially designed for Fire and Rescue services. It can hold any PAC ventilator in a special compartment, cylinder, regulator, hoses and masks, with plenty of space for First Aid kit, stiff neck collars and any other medical equipment. The bag is flashed with a recognisable Trauma logo on both ends. It has a replaceable base and can be stand horizontal or vertical for easy storage.



Rescue Pack Storage Bag

The Rescue Pack bag is specially suitable for Fire and Rescue services. It has multiple storage facilities to hold a PAC ventilator, cylinder, regulator, hoses and first flow masks, with specially designed side and front pockets for other medical equipment.

Specification

Material: Heavy duty Cordura, washable and water resistant

Labels: Removeable side labels with Velcro fastenings and clear pockets to slide in labelling. Flashed with Pneupac logo as standard. (Company) logo can be configured to customer requirements for which an additional cost may be charged.

Zips: Y2K high quality heavy duty zips

Dimensions (approx):

Trauma Pack Bag: 600mm W x 280mm H x 300mm D

Rescue Pack Bag: 560mm W x 300mm H x 320mm D

Ordering Information

Medical Storage Bags	Order No
Trauma Pack Bag	W1236
Rescue Pack Bag	W1226

Regulators

Ordering Information

Oxygen	Order No
Oxygen Pin-Indexed Cylinder Connection Regulator with single Outlet for BS5682 Probe	500-A162/CE
Oxygen Pin-Indexed Cylinder Connection Regulator with single Outlet for BS5682 Probe & Variflow Selector Settings: (0, 5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10 & 15 LPM) <i>Note Vari-flow Selector Valve Accessories A928/5 - Oxygen Tube & A369 Therapy Mask will also be required</i>	500-A162CE/Z43
Oxygen Pin-Indexed Cylinder Connection Regulator with Variflow 0- 15 LPM for Oxygen Therapy only	500-A162/Z20/CE
Bullnose Cylinder Connection Regulator with single Outlet for BS5682 Probe	500-A166/Z31/CE

Air	Order No
Air Pin-Indexed Cylinder Connection Regulator with single Outlet for BS5682 Probe	500-A162CE/228

Entonox*	Order No
Entonox Pin-Indexed Cylinder Connection Regulator with single Outlet for BS5682 Probe	500-A162/Z5/CE
Entonox Pin-Indexed Cylinder Connection Regulator with 1/4 BSP Outlet for Inhalation Analgesia Sets <i>* Entonox is a BOC Trademark</i>	500-A162/Z6/CE

Accessories

	Order No
Gauge, Contents for Entonox Regulators	W6829/CE/Z1
Oxygen Therapy Mask, Disposable (pack of 5)	500-A369/5
Oxygen Therapy Tube for use with Disposable Therapy Mask (pack of 5)	500-A928/5
Shroud, Gauge, Shock Absorbing	500-864
Washers, Pin Index Yoke (Pack of 3)	500-276
Wheel, Valve, Cylinder	500-A171

Ordering Example: A162/CE/Z43 + 500-A369/5 + 500-A928/5

- = **1)** Pin index regulator with single outlet & variflow (500-A162CE/Z43)
- + **2)** Disposable Oxygen therapy masks (500-A369/5)
- + **3)** Oxygen therapy tubes for use with 500-A369/5 (500-A928/5)